



»Lead the way to
quick and clear H&E staining«

MMI H&E STAINING KIT PLUS

The only H&E staining kit designed for contamination free staining.
For users with small and medium number of stained slides a day!



Content:

15 mmi SafeStain Ampoules Harris Hematoxylin

15 mmi SafeStain Ampoules Eosin aqueous

1 manual

Customer benefits:

Absolute RNase free

Eliminate contamination by using the mmi SafeStain Ampoules.

Quick and safe

Reduce RNA digestion while staining. Use the quick 4 min mmi Staining Protocol.

All tissues

Applicable for fresh frozen and FFPE samples.

Small and ecological


Single use mmi SafeStain Ampoules. Eliminate wasting H&E. Just use one mmi SafeStain Ampoule to stain 3-4 slides.

Smart and economical

MMI SafeStain Ampoules are cost effective for users staining less than 30 slides a day.

Clean and Easy

The handy mmi SafeStain Ampoules help you to stain your samples without the need of stain filled jars.



Best quality combined
with high safety and
comfort!



MMI AG
Flughofstrasse 37
8152 Glattpfurgg, Switzerland
Phone: +41 44 809 10 10
Fax: +41 44 809 10 11
Mail: info@molecular-machines.com
www.molecular-machines.com

MMI GmbH
Breslauer Strasse 2
85386 Eching, Germany
Phone: +49 89 319 048 40
Fax: +49 89 319 048 59
Mail: info@molecular-machines.com
www.molecular-machines.com

MMI Inc.
P.O. Box 348
Haslett, MI 48840, USA
Phone: +1 603 629 9536
Fax: +1 321 978 0304
Mail: sales_us@molecular-machines.com
www.molecular-machines.com